

Pose Estimation in NWB

ndx-pose



Neurodata
Extensions
Catalog



Pose Estimation in NWB

ndx-pose

NWB Currently supports writing from two popular pose estimation frameworks

- SLEAP
- DeepLabCut (DLC)

We do not have time to cover exactly *how* to use or understand these packages

- See the main keynote talk by Mackenzie Mathis instead

So we will simply explain the input-output (I/O) workflow

Note: Storing natural behavior in NWB

- The NWB data type for videos called the *ImageSeries*
- It has two modes
 - Internal
 - each frame of the video is extracted and written as a Dataset
 - same as any other array-valued series you might find in an NWB file
 - External
 - data type instead contains a *path* to the video file(s) outside the NWB file
- Reasoning for this is a bit nuanced...
- For history of decision, see...
 - [NWB:#1647](#), [DANDI:#769](#) and [DANDI-helpdesk:#30](#)

Note: Storing natural behavior in NWB

■ Best Practices

- [Internal mode](#)
 - videos of physiological recordings (such as optical imaging)
 - stimulus presentations
- [External mode](#)
 - videos of natural behavior

■ DANDI

- strongly recommends these practices for storage
- will automatically remap all paths and organize folder storage at time of upload

Step 1: Install SLEAP or DLC

- SLEAP - Talmo Lab - [instructions](#) - [GitHub repo](#)
- DLC - Mathis Lab - [instructions](#) - [GitHub repo](#)

- This step can take a while, but just be patient and persistent

- If anything goes differently from instructions, when in doubt ask
 - Google
 - ChatGPT
 - or raise an issue on their repos/helpdesks

Step 2: Launch GUI

- These can be run headless through their API but I recommend using the **G**raphical **U**ser Interface (GUI) to help walk you through it
- SLEAP can import NWB ImageSeries and export NWB PoseEstimation
- Both are able to load in original videos (.mpg, .mp4, .avi, etc...)
- Due to the previous discussion of storing ImageSeries of natural behavior it is recommended to always load in the raw video to both platforms

Steps 3-7

- Step 3: Load videos
- Step 4: Extract frames for training
- Step 5: Apply points and labels on training frames
- Step 6: Train the model
- Step 7: Run the prediction

Saving results

- The output of the prediction over the entire video can then be saved in either the source `.slp` or `.h5` (DLC) formats
- These can in turn be converted to NWB format via two utility packages
 - `Dlc2nwb` (`.h5` \rightarrow `.nwb`)
 - `sleap-io.nwb` sub-module (`.slp` \rightarrow `.nwb`)
- BONUS: SLEAP can export directly to NWB through the GUI

Understanding PoseEstimation

- In either case (or for any future integrations with keypoint tracking tools) these are modelled using a common structure in NWB called the `PoseEstimation`
- A single `PoseEstimation` is a container of multiple `PoseEstimationSeries`, one for each node (usually a body part) being tracked
- The `data` field of each `PoseEstimationSeries` represents the `x/y/z` positions being tracked over time
- There is also the `confidence` ('likelihood') associated with each position over time

Understanding PoseEstimation

- See it used in practice in the IBL Brain Wide Map (dandiset #409)

```
nwbfile.processing["behavior"]["PoseEstimationLeftCamera"]
```

PoseEstimationLeftCamera abc.PoseEstimation at 0x140450551033632

Fields:

description: Estimated positions of body parts using DeepLabCut.

nodes: <StrDataset for HDF5 dataset "nodes": shape (11,), type "<O">

pose_estimation_series: {

 nose_tip <class 'abc.PoseEstimationSeries'>,
 paw_l <class 'abc.PoseEstimationSeries'>,
 paw_r <class 'abc.PoseEstimationSeries'>,
 pupil_bottom_r <class 'abc.PoseEstimationSeries'>,
 pupil_left_r <class 'abc.PoseEstimationSeries'>,
 pupil_right_r <class 'abc.PoseEstimationSeries'>,
 pupil_top_r <class 'abc.PoseEstimationSeries'>,
 tongue_end_l <class 'abc.PoseEstimationSeries'>,
 tongue_end_r <class 'abc.PoseEstimationSeries'>,
 tube_bottom <class 'abc.PoseEstimationSeries'>,
 tube_top <class 'abc.PoseEstimationSeries'>

}

source_software: DeepLabCut

```
nwbfile.processing["behavior"]["PoseEstimationLeftCamera"]["nose_tip"]
```

nose_tip abc.PoseEstimationSeries at 0x140450549420416

Fields:

confidence: <HDF5 dataset "confidence": shape (242446,), type "<f8">

conversion: 1.0

data: <HDF5 dataset "data": shape (242446, 2), type "<f8">

interval: 1

offset: 0.0

reference_frame: (0,0) corresponds to the upper left corner when using width by height convention.

resolution: -1.0

timestamps: <HDF5 dataset "timestamps": shape (242446,), type "<f8">

timestamps_unit: seconds

unit: px